



## Production of X-Rays



## History



**Wilhelm  
Conrad  
Roentgen**

**Discovered  
x-rays on  
Nov. 8, 1895**



## Dr. Otto Walkhoff

**First Dental Radiograph  
25-minute exposure  
( $< 0.5$  seconds today)**




before exposure



after exposure


## Dr. C. Edmund Kells

**First Intraoral Radiograph  
1896**




# X-ray Physics

## Atom (electrically neutral)



Atomic Number (Z) = # of protons




## ELECTROSTATIC FORCE

Attraction between protons and electrons



## CENTRIFUGAL FORCE

Pulls electrons away from nucleus



Balance between electrostatic force and centrifugal force keeps electrons in orbit around nucleus

## Binding Energy

The amount of energy required to remove an electron from its orbit (= electrostatic force). Depends on atomic number (# of protons).


## Electromagnetic Radiation

The movement of energy through space as a combination of electric and magnetic fields

Travel at the speed of light  
( $3 \times 10^8$  meters/second)  
(186,000 miles/second)


## Electromagnetic Radiation

- X-ray
- Visible light
- Gamma ray
- TV waves
- Microwaves
- Radio waves



Wavelength x Frequency = Speed of light

## Electromagnetic Spectrum




## ENERGY

Ability to penetrate

Shorter wavelength, higher energy

Higher frequency, higher energy



## X-ray Characteristics


- High energy waves
- No mass
- No charge (neutral)
- Travel at speed of light
- Invisible

## X-ray Characteristics


- Travel in straight line
- Cannot be focused to a point
- Differentially absorbed
- Cause fluorescence
- Harmful to living tissue






## X-ray Machine

**Tubehead**  
**Support arms**  
**Control panel**




PID = position indicating device  
BID = beam indicating device






## Thermionic Emission

Release of electrons  
from hot filament




Increase current, increase heat  
Increase heat, increase # electrons




## Anode

Dental x-ray machines have stationary anode





## Medical Radiography



## Tungsten (Filament and Target)

- High atomic number ( $Z=74$ )
- Transfers heat readily
- High melting point ( $3370^{\circ}\text{C}$ )
- Can be drawn into fine wire


## X-ray Tube Components




The smaller the focal spot (target), the sharper the image (teeth) will be.

During x-ray production, a lot of heat is generated. If the target is too small, it will overheat and burn up.


## Line Focus Principle




## X-ray Machine Components




## 60-cycle AC



## 60-cycle Alternating Current





$$\frac{\text{Number of Impulses}}{60} = \text{Seconds}$$

$$60 \text{ impulses}/60 = 1 \text{ second}$$

$$30 \text{ impulses}/60 = 0.5 \text{ second}$$

$$15 \text{ impulses}/60 = 0.25 \text{ second}$$




## Exposure Switch



Allows current to flow to complete high and low voltage circuits.

You cannot overexpose by holding the exposure switch down too long!




milliAmpere (mA) selector

## Step-Down Transformer

220 volt → 3 – 5 volts

Filament Circuit

## Step-Down Transformer






## Autotransformer


Controls voltage between  
anode and cathode

Regulated by kVp selector

(Similar to a rheostat)



## Autotransformer





## Step-Up Transformer


65 - 90 V  $\rightarrow$  65,000 - 90,000 V

(65 kVp - 90 kVp)

kVp = kiloVoltage peak

## Step-Up Transformer






## X-ray Production

### Bremmstrahlung

### Characteristic




### Bremsstrahlung Radiation


(Braking radiation, general radiation)

X-rays produced when high-speed electrons from the filament are slowed down as they pass close to, or strike, the nuclei of the target atoms

### Bremsstrahlung X-ray Production



### Bremsstrahlung X-ray Production Maximum energy



### Characteristic Radiation

X-rays have energies characteristic of the target material (energy = difference between binding energies of target electrons involved, e.g., K & L, K & M, etc.)

The energy of the high-speed electron from the filament must be higher than the binding energy of the target electron with which it interacts in order to eject the target electron


## Tungsten

K-shell binding energy = 70 keV



L-shell binding energy = 11 keV

M-shell binding energy = 3 keV


## Characteristic X-ray Production



## Characteristic X-ray Production




What is the energy of the characteristic x-ray produced? (M-shell binding energy = 3 keV)



What is the energy of the characteristic x-ray produced above?

## X-ray Spectrum (variable x-ray energies)




## X-ray Spectrum results from:



Varying electron/nucleus distances



Multiple electron interactions



Varying voltage (AC)




Electron interactions with target atoms

1 % produce x-rays

99 % produce heat

The excess heat is controlled by high melting point of tungsten, conductive properties of copper sleeve, and cooling from oil surrounding x-ray tube

## X-ray Beam Modifiers




## Exposure Factors


kVp

mA

Exposure time



Incorrect exposure factors  
(too dark)




Correct exposure factors



Incorrect exposure factors  
(too light)


kVp (kiloVolt peak)




**Increasing kVp results in:**

- Higher average energy of x-rays
- Greater maximum energy x-rays
- More x-rays

**mA (milliampere)**  
(kVp, exposure time constant)



**Exposure time**  
(kVp, mA constant)



**Increasing mA or Exposure Time results in:**

**An increase in the number of x-rays produced**

**No change in the energy of the x-ray beam**

**mAs or mAi**

**milliamperes (mA) x seconds (s)**

**milliamperes (mA) x impulses (i)**

**60 impulses = 1 second**


**10 mA x .5 seconds = 5 mAs**

**20 mA x .25 seconds = 5 mAs**

**mAi = 60 x mAs**

**Constant patient size**


1. Proper kVp, mA, exposure time (e.t.) **B**
2. Increase mA; no change in kVp, e.t. **A**
3. Decrease e.t.; no change in kVp, mA **C**
4. Increase kVp; no change in mA, e.t. **A**
5. Double mA, halve e.t.; no change in kVp **B**



## Filtration

The process of removing low-energy x-rays from the x-ray beam

## Total Filtration





The filter is usually located in the end of the PID which attaches to the tubehead.




## Collimation

Regulates the size and/or shape of the x-ray beam

- ↓ area covered (less patient exposure)
- ↓ scatter radiation





**Quality**

average energy

**Quantity**

number of x-rays

	Quality	vs.	Quantity
kVp	↑ (1°)		↑ (2°)
mA	No change		↑
Time	No change		↑
Filtration	↑		↓

**Half-Value Layer**

Indicates the quality (energy) of the x-ray beam


The thickness of aluminum needed to reduce the number of x-rays by one-half


2.5 - 70 kVp+; 1.5 < 70 kVp

**INVERSE SQUARE LAW**

The intensity of radiation varies inversely as the square of the target \* -film distance


\* target = source, focal spot, focus






**D = distance**  
**KI = known intensity**  
**UI = unknown intensity**

## Inverse Square Law



1. Divide distance 1 (known intensity) by distance 2 (unknown intensity)
  2. Square this number
  3. Multiply this number by the known intensity
- This will equal the intensity at distance II


If the intensity of an x-ray beam is 4 at 12 feet, what is the intensity of the x-ray beam at 3 feet?



(I) at 12 ft. = 4    (I) at 3 ft. = ?

$12 \div 3 = 4$   
 $4 \text{ squared} = 16$   
 $16 \times 4 = 64$  (Intensity at 3 feet)

If the intensity of an x-ray beam is 72 at 2 feet, what is the intensity of the x-ray beam at 12 feet?





**(I) at 2 ft. =72    (I) at 12 ft. = ?**

$$2 \div 12 = 1/6$$

$$\left( \frac{D_{KI}}{D_U} \right)^2 \times KI = UI$$

$$1/6 \text{ squared} = 1/36$$

$$1/36 \times 72 = 2 \text{ (Intensity at 12 feet)}$$