

## Pontics and the Edentulous Ridge

Donna N. Deines, DDS, MS

Resources: Shillingburg et al: *Fundamentals of Fixed Prosthodontics*  
 Rosenstiel et al, *Contemporary Fixed Prosthodontics*  
 Eissmann, H.F. in *Dental Laboratory Procedures: Fixed Partial Dentures* (Rudd, Morrow)  
 Chiche, G & Pinault, A: *Esthetics of Anterior Fixed Prosthodontics*

## Pontic Design and the Edentulous Ridge



### Esthetic

Appearance of replacement  
 "Emergence" from ridge  
 Space for porcelain

### Biologic


Ease of cleaning  
 Patient comfort  
 Healthy tissue  
 Occlusion

### Mechanical

Rigid: resist deformation /  
 porcelain fracture  
 Strong connectors




## Pontic Materials and Design



- Glazed porcelain
- Highly polished gold
- M-C finish line: not on edentulous ridge
- Acrylic / composite resin (porous) – poor surface
- Proper design most important




## Alveolar Ridge Resorption



Normal Minimal Excessive



Edentulous space defined by its boundaries:

- Edentulous ridge
- Opposing occlusion
- Tongue and cheek





Eissmann, Harold F.

## Alveolar Ridge Resorption (Minimal)




- Broad ridge / vertical height
  - Uncommon
  - Can compromise occlusocervical dimension →
    - decreased strength of restoration
    - limits restoration possibilities

## Changes in the edentulous ridge




- Alveolar resorption and remodeling
  - Apical and facial resorption
  - Greater with trauma / periodontal disease
- Pontic or tissue modifications for esthetics / cleansability



### Residual Ridge Deformities

Siebert, JS (*Compend Contin Educ Dent* 4:437, 1983)




Normal: (8%)      Class I: F-L Width (32%)

Class II: O-G Height (3%)      Class III: F-L and OG Height (56%)

Lateral resorption → narrow ridge  
Vertical resorption → ridge height defect


Shillingburg

### Pontic Contour Changes with Ridge Resorption



- A pontic should have the same incisogingival height as original (or contralateral) tooth.
- With resorption, the pontic becomes longer in order to contact the ridge concavity.
- *Contour must be blended smoothly to avoid a ledge at the cervical.*

### Ridge Resorption – Elongated Pontic




### Esthetic Modification: Gingival Contour of Pontic




- Contour in apical ½ should approximate length of adjacent teeth.

### Esthetic Modification: Gingival Contour of Pontic




- Contour in apical ½ should approximate length of adjacent teeth.
- The facial surface must be altered to curve gently from the G-F angle to the middle of facial surface

### Esthetic Modification of Pontic Contour



### Esthetic Modification of Pontic Contour




### Biological Considerations: Tissue Contact



- Any contact should be pressure-free
  - Tissue should not blanch
- Keratinized attached tissue
- No contact with ridge if possible or smallest possible area with convex pontic surface

### Tissue Contact of Pontic




Contact on non-keratinized mucosa / low frenum attachment must be adjusted to avoid ulceration.

### Ridge-Lap and Tissue Impingement



- Pressure by pontic → tissue necrosis
- Keratinized tissue contact only
- Avoid non-keratinized tissue and frenum attachments.
- Floss should pass through without resistance.

### Tissue Contact of Pontic (Ridge Lap Pontic)




- Severe ridge lap of pontic prevents cleaning.
- Pontic pressure causes tissue ischemia and necrosis.

### Tissue Contact of Esthetic Pontics





- Area of contact small and convex.
- No space between pontic and soft tissue of facial side of ridge
- Contact on attached keratinized gingiva only
- No pressure on the ridge

## Oral Hygiene




- Convex smooth surface of pontic
- Adequate gingival embrasures
- Patient education and hygiene aids
  - Floss threader; Super-floss / gauze; Proxy-brush; WaterPik



- Connectors and pontic contours must not impinge on gingival tissue → gingival inflammation / recession / necrosis.
- Yet there must be adequate thickness for strength
- Strength properties of materials
  - Metal / porcelain / acrylic




## Mechanical Considerations




- Strength: O-G thickness of pontic / connector
- Material: metal > metal-ceramic > ceramic
- Position: straight line as possible
- Occlusion: normal centric contacts
  - Occlusal table commensurate with occlusion
  - Narrowing pontic ≠ decreased occlusal force

## Surgical Ridge Correction:




**Cuff of tissue** adjacent to edentulous space  
Removal allows ↑connector size / cleaning



## Pontic Designs




## Hygienic Pontic (Non-Esthetic)




- (PRO) Easily cleaned
- (CON) Entrapment of food debris
- No contact with the residual ridge
- ~2-3mm space between apical portion of pontic and ridge
  - (tissue proliferation if too close)
- At least 3mm thick for strength (O-C)
- Convex contour FL / MD ("fishbelly")

### Modified Hygienic Pontic ("Perel")


- Increased strength in connectors
- Decreased deflection with decreased O-C height



### Conical Pontic




- Small convex area of contact at center of ridge.
  - Easy to clean / comfortable for patient
- Facial / lingual contours depend on residual ridge width




**Contraindication:**  
Point contact w/ broad ridge → food entrapment

### Esthetic Pontic Designs


- Ridge Lap
- Modified Ridge Lap
- Ovate



### Modified Ridge Lap Pontic




**Ridge-lap Pontic:**  
Not cleansable




### Modified Ridge Lap Pontic




- Ridge contact on the facial half of ridge
- Contour of tissue-contacting area is convex
- (Slight faciolingual concavity on facial side of ridge)
- Tissue contact resembles a "T"




Intimate tissue contact on facial half helps prevent debris accumulation



### Tissue Contact of Maxillary FPD




### Ovate Pontic Design




- Blunt / rounded apex set into ridge concavity
- Requires broad flat (rounded) ridge
- Concavity formed by:
  - Immediate provisional restoration
  - Surgical creation


### Surgical Correction for Ovate Pontics: Gingivoplasty - ridge concavity / pontic convexity



All-Ceramic Resin-Bonded FPDs  
#6-8 and #9-11  
Lawrence Brecht, DDS  
ACP Prosthodontia



### Ovate Pontic: Immediate Provisional FPD



- Facial contour extends 2.5 mm apical to FGM
  - reduced to 1mm – 1.5mm at 4 weeks

### Ovate Pontic: Immediate Provisional FPD





May wait ~6 months before tissue stabilizes and final restoration is fabricated.



### Pre-Treatment Assessment


- Pontic Space
  - M-D space (tipping)
  - I-G space (supraeruption)
  - Orthodontics
  - Diagnostic Waxing
- Residual Ridge Contour
  - Pontic modifications
  - Surgical modifications
  - Gingival architecture preservation





Form is compromised in the lesser visible half.


### Awkward Pontic Spaces



- Duplicate the dimension of the more visible mesial half of the adjacent tooth.
- (B) Narrow space
- (C) Wide space


*(Redrawn from Blancheri RL: Rev Assoc Dent Mex 8:103, 1950.) in Rosenstiel: Contemporary Fixed Prosthodontics.*

### Pontic Modifications: "Black Triangles" – open embrasures / resorption



- Unesthetic loss of interdental papillae with ridge resorption
- Plaque accumulation / phonetic difficulties
- Interference with flossing

### Pontic Modifications: Excessive Length



- Vertical resorption requires excessive pontic length
- Contour retainers / pontics normally; shape and stain pontic extension to simulate root surface.

### Esthetic Modifications of Pontic Contour



- Apparent shortening of crown with cervical shading
- Pink porcelain to simulate gingivae
- Decreased embrasure space with cervical contouring

### Pontic Modification: Embrasure filled with porcelain Narrow faciolingual ridge width - Mandibular




- Esthetic problems (?)
- Cleansability usually not a problem


### Soft Tissue Augmentation Surgery



Hard tissue grafts - implants  
Connective tissue grafts:  
(subepithelial / submucosal)




### Soft-Tissue Augmentation Surgery



Vertical resorption would require excessive pontic length of centrals (length : width too great)

Chiche & Pinault

### (Post-surgical) Soft Tissue Augmentation Restoration




Balance

- Tooth length:width
- Central:Lateral proportion

Chiche & Pinault

### Management of deficient ridge with shaded gingival porcelain




Malament, KA and Neeser, S  
Dent Clin N Am, 2004 735-44


FPD #6-11 inserted 1982

Residual ridge resorption presented 2000

### Pink Porcelain Simulation Gingiva (Implant Abutments)



### "Andrews Bridge"



Fixed retainers are connected by a rectangular bar which follows the curve of the ridge.

Denture teeth in an acrylic resin flange that clips over and is stabilized by the bar.