

Paralleling Technique




Patient Preparation:

Seat patient

Adjust headrest

Place lead apron, thyroid collar



Long axis of the tooth parallel with long axis of the film

Paralleling Technique

ADVANTAGES

(compared to Bisecting Angle technique)


1. Better dimensional accuracy
2. Beam alignment simplified
3. Easier to standardize films
4. Head position not critical

Paralleling Technique







DISADVANTAGES







(compared to Bisecting Angle technique)

1. Uncomfortable
2. Limited by anatomy




Film/tooth/ring all parallel
Central ray perpendicular to tooth/film







Posterior Periapical




The film should be equidistant from the teeth in an anterior-posterior direction




The teeth being radiographed must be in contact with biteblock.








Cotton rolls may be used in any area of the mouth to help support the biteblock, especially if an edentulous region opposes the teeth being radiographed; often makes it more comfortable for the patient. The cotton roll should be placed against the arch opposite to the one being radiographed.





While supporting the bar with the fingers of one hand, slide the ring down close to the face with the other hand.




Make sure head is supported by headrest before exposing films.




Maxillary Canine



Note overlap between canine and first premolar. Usually not avoidable in maxillary canine region using paralleling technique.




Maxillary Premolar




Film positioned so that anterior edge is at least in the middle of the canine, or more anterior if possible. Approximately centered on 2nd premolar. Top edge of film in center of palate.

Maxillary Premolar





Maxillary Molar




Film centered on second molar. Top edge of film in center of palate.


Maxillary Molar



Mandibular incisor





Film centered on midline. Film should be placed away from the teeth as much as possible.




Let patient guide film into place as they close their mouth

Mandibular incisor







Mandibular Canine




Film centered on canine. Film should be placed away from the teeth as much as possible.

Mandibular Canine




Mandibular Premolar





Anterior edge of film is positioned at least in the middle of the canine, or more anterior if possible. Approximately centered on 2nd premolar. The film should be placed more toward the middle of the mouth, away from the teeth. This is usually the most uncomfortable film taken on a patient using the paralleling technique.

Mandibular Premolar






Mandibular Molar



Film centered on second molar. Film may be placed next to the teeth.


Mandibular Molar




For patients with long roots:

Normal alignment
misses apices

Increase vertical
Raise PID










top edge of PID
above ring








Paralleling Technique

ERRORS









Paralleling Film Placement

- Film size: # 1 ant., # 2 post.
- Film holder: Rinn XCP
- White-on-white
- Dot-in-the-slot
- Vert. anterior, hor. posterior
- Correct ant.-post. placement

Paralleling Film Placement

- Equal tooth-film distance
- Biteblock in contact
- Cotton roll for stability (opposite arch)
- Ring close to skin